

ISPGR WS



Bridging the Gap

Best Practices in Mobile Brain Imaging

Pre-processing EEG gait data



Differences between biomechanical and EEG data

- Heel and toe marker kinematics (low-pass filtered)






Differences between biomechanical and EEG data

- Heel and toe marker kinematics (low-pass filtered)






Differences between biomechanical and EEG data

- Heel and toe marker kinematics (low-pass filtered)
 - Raw data



Differences between biomechanical and EEG data

- Heel and toe marker kinematics (low-pass filtered)
 - Filtered data



Objectives of this session

After this session you will

- know the steps involved in EEG preprocessing
- know the different type of artefacts that occur during (mobile) EEG and,
- know how to determine the physiological origin of artefacts and brain data

What is preprocessing?

- All the steps required to convert the raw data to interpretable data
- Iterative process consisting of (at least) the following steps:
 - Visual inspection of the data
 - Band pass filtering
 - Removal of noisy data (and re-referencing)

What is preprocessing?


- All the steps required to convert the raw data to interpretable data
- Iterative process consisting of (at least) the following steps:
 - Visual inspection of the data
 - Band pass filtering
 - Removal of noisy data (and re-referencing)
 - **Independent component analysis for artefact identification**

What is preprocessing?

- All the steps required to convert the raw data to interpretable data
- Iterative process consisting of (at least) the following steps:
 - Visual inspection of the data
 - Band pass filtering
 - Removal of noisy data (and re-referencing)
 - **Independent component analysis for artefact identification**
- Open-source toolboxes for Matlab (fieldtrip, eeglab) or Python (MNE-python)


Why preprocessing?

- Unavoidable artefacts during mobile brain imaging
 - Scalp and neck muscle activations
 - Eye movement and blinks
 - ECG artefact
 - Movement artefacts (i.e., time locked to gait events such as heel strikes)
- Environmental noise (i.e., line noise at 50 or 60 Hz and harmonics, WiFi, bluetooth)



Why preprocessing?

- Unavoidable artefacts during mobile brain imaging
 - Scalp and neck muscle activations
 - Eye movement and blinks
 - ECG artefact
 - Movement artefacts (i.e., time locked to gait events such as heel strikes)
- Environmental noise (i.e., line noise at 50 or 60 Hz and harmonics, WiFi, bluetooth)



Step 1: visual inspection


- First step: visually inspection of the data to:
 - Get an impression about the data quality
 - Verify that all events/conditions are correctly specified in the recording
 - Remove irrelevant data segments (i.e., breaks between consecutive experimental tasks)

Step 1: visual inspection

- First step: visually inspection of the data to:
 - Get an impression about the data quality
 - Verify that all events/conditions are correctly specified in the recording
 - Remove irrelevant data segments (i.e., breaks between consecutive experimental tasks)
- Preprocessing will not be able to fix low quality data or recover relevant event information absent in the raw data


Step 2. Identification and removal of bad channels

- Visual inspection on data quality:
 - Are there any channels that are noisy/flat throughout the recording?





Step 2. Identification and removal of bad channels

- Visual inspection on data quality:
 - Are there any channels that are noisy/flat throughout the recording?
- If a channel is noisy/flat for the entire recording -> remove it
- If noise spreads over multiple channels for short time windows > remove the segment





Step 3. Filtering

- Why filter EEG data?
 - Remove low-frequency drift
 - Remove high-frequency noise





Step 3. Filtering

- Why filter EEG data?
 - Remove low-frequency drift
 - Remove high-frequency noise



Step 3. Filtering

- Why filter EEG data?
 - Remove low-frequency drift
 - Remove high-frequency noise




4. Re-referencing

- EEG is recorded relative to a reference point
- The ‘common average’ reference is most frequently used in high density EEG recordings:

$$\bullet V_i = V_i - \frac{1}{N} \sum_{j=1}^N V_j$$


Average over all N
EEG channels

- If a channel is removed from the data, re-reference the data



5. Remove noisy data segments


- Remove noisy segments of data
 - Noisy data during the experiment (muscle bursts/movement artefacts)
 - Prolonged periods of recorded EEG during breaks between experimental tasks
- Create short epochs of 1 second to avoid the unnecessary removal of clean data



6. Independent component analysis

- A statistical technique to separate multivariate mixture signals into **maximally independent** components(ICs) using higher-order statistics
- Solving the cocktail party problem


The Cocktail Party Problem
SOLVING WITH ICA



6. Independent component analysis

- A statistical technique to separate multivariate mixture signals into **maximally independent** components(ICs) using higher-order statistics
- Solving the cocktail party problem
- **The brain as a cocktail party**
- Each electrode measures a linear combination of brain sources


The Cocktail Party Problem
SOLVING WITH ICA



6. Independent component analysis


- A statistical technique to separate multivariate mixture signals into **maximally independent** components(ICs) using higher-order statistics
- Solving the cocktail party problem
- **The brain as a cocktail party**
- Each electrode measures a linear combination of brain sources
- The mixing matrix describes how an ICs project to the scalp electrodes

The Cocktail Party Problem **SOLVING WITH ICA**



Independent component characteristics

Electrode locations




Independent component characteristics

Electrode locations




Projection of a single component




Independent component characteristics


Electrode locations



Projection of a single component




Modelling component as a single dipolar source





Independent component characteristics


Electrode locations




Projection of a single component



Modelling component as a single dipolar source




Activity power spectrum




Independent component characteristics


Electrode locations



Projection of a single component





Modelling component as a single dipolar source





Artefact classification

- Based on:
 - Scalp map: describes how an independent component distributes to individual channels
 - Time- and frequency domain representation
 - Dipole location (inside or outside the head) and goodness of fit



Independent component classification (1)





Independent component classification (1)



Independent component (2)



Independent component (2)




Independent component classification

- Scalp maps
- Time domain activity
- Power spectral density
(frequency domain)
- Dipole location and goodness
of fit


Independent component classification

- Scalp maps
- Time domain activity
- Power spectral density (frequency domain)
- Dipole location and goodness of fit
- **Classify all components**




Independent component classification

- Scalp maps
- Time domain activity
- Power spectral density (frequency domain)
- Dipole location and goodness of fit
- **Classify all components**
- **ICLabel plugin**



Artefact classification

- Scalp maps
- Time domain activity
- Power spectral density



What about movement artifacts?

of fit


- Classify all components
- ICLLabel plugin

Movement artefact classification

- 1) Define trials around a specific gait event of interest, i.e., heel strike
- 2) Calculate the consistency of independent components over trials (i.e., SNR)


Movement artefact classification

- 1) Define trials around a specific gait event of interest, i.e., heel strike
- 2) Calculate the consistency of independent components over trials (i.e., SNR)




Artefact classification


- Scalp maps
- Time domain activity
- Power spectral density (frequency domain)
- Dipole location and goodness of fit
- **(gait) event-related potential**



Preprocessing prepares data for analysis



Preprocessing prepares data for analysis



Next step: Processing of EEG data in relation
to the gait cycle

Thank you for your attention!