



ISPGR WS
Bridging the Gap
Best Practices in Mobile Brain Imaging
Maastricht 2025

EEG signal




Traditional EEG experiments




MOBI EEG experiments

Option 1: Lab task during walking




Option 2: Lab task during walking



**How can I record EEG, EMG,
MoCap, video, ... in a
synchronized manner?**

Part 1- How to synchronize recordings



LabStreamingLayer (LSL)

LSL distribution comprises:

- **Core Library:** liblsl (C, C++, Python, Java, MATLAB).
- **Platform:** Cross-platform (Windows, Linux, macOS, Android).

The most common way to use LSL is to use one or more applications to stream data from one or more devices (e.g., EEG and EMG) over the local network and record the with the LabRecorder.

Steps to use LSL:

1. Per device create LSL outlet
2. Fetch data from device
3. Push data to LSL


Example Matlab Code

```
1 %% instantiate the library
2 disp('Loading library...');
3 lib = lsl_loadlib();
4
5 % make a new stream outlet
6 disp('Creating a new streaminfo...');
7 info = lsl_streaminfo(lib,'DelSys','EMG',8,100,'cf_float32'); % 8 channels, 100 Hz, float32
8
9 % initiate DelSys Trigno wireless EMG SDK
10
11 DelsysInput = tcpip(HOST_IP,50041); % HOST_IP is the IP address of the computer running the Delsys SDK, 50041
12 DelsysInput.InputBufferSize = 6400; % Buffer size for the input stream
13
14 disp('Opening an outlet...');
15 outlet = lsl_outlet(info);
16
17 % send data into the outlet, sample by sample
18 disp('Now transmitting data...');
19 while true
20     % get data from device
21     tmp_data = fread(DelsysInput,bytesReady); % read data from device (8 channel EMG)
22
23     % push data to LSL outlet
24     outlet.push_sample(tmp_data); % push data to LSL outlet
25 end
```




EEG LSL outlet

`lsl_streaminfo(lib,'BioSemi','EEG',8,100,'cf_float32')`



EMG LSL outlet

`lsl_streaminfo(lib,'TMSi','EMG',12,200,'cf_float32')`



LSL lib


EEG Stream



EMG Stream



EMG




Buffer

pulls data



EMG LSL outlet

writes data



00.10291 s

00.10450 s

00.10670 s


..

..

..

.xdf

Why is subsecond precision important?



Take home

1. EEG plus other modalities are very interesting!
2. LSL can be used to record multimodal data in a highly synchronized manner
3. Sub-second precision is important for MOBI with EEG!

Thank you for listening 😊