

fNIRS and Gait

Sue Peters PT, PhD

Assistant Professor, School of Physical Therapy, UWO

Scientist, Gray Centre for Mobility and Activity, Lawson Health Research Institute

Associate Member, Center for Brain and Mind

Registered Physiotherapist



Western
HealthSciences






Western
The Brain and
Mind Institute



THE GRAY CENTRE FOR
MOBILITY AND ACTIVITY
at Parkwood Institute

Brain Anatomy




Human Brains at Rest and During a Task

At Rest

- Brain is 2% of body weight, but receives 20% of cardiac output
- Brain blood flow/metabolism are generally constant

During a task

- Increased oxygen consumption and brain blood flow
- Importantly → Increases are brain region specific



(Review: Lauritzen et al Neuroimage. 2012)


Human Brains at Rest and During a Task

Proposed Mechanisms

- Neurovascular coupling
 - Oxygen metabolism to support neuronal activity
- Overall balance toward excitation or inhibition
 - Generates increase in regional brain blood flow


Hemodynamic response

- Oxy-Hemoglobin (HbO)
- Deoxy-Hemoglobin (HbR)
- Figure → typical to see HbO in red and HbR in blue
 - Here include 2 standard error around the mean




functional Near Infrared Spectroscopy (fNIRS)


A



Review (Chen et al. Front Neurosci 2020)



functional Near Infrared Spectroscopy (fNIRS)



Review (Chen et al. Front Neurosci 2020)



NIRSport 2 system we will use today

Advantages and Disadvantages of fNIRS

- Advantages
 - Resistant to motion artifact
 - Portable (can be done at bedside, outdoors, etc)
 - Quick set up (~5 minutes) and clean up
 - Low cost to operate (vs MRI)
 - Better temporal resolution than fMRI
 - fNIRS measures oxy- *AND* deoxy- Hb (and total Hb)
 - Better spatial resolution than EEG
 - Easily combined with other modalities
 - **Field is changing very rapidly**
- Disadvantages
 - Worse temporal resolution than EEG
 - Worse spatial resolution than fMRI
 - **Field is changing very rapidly**



NIRSport 2 system we will
use today

EEG vs fNIRS

- Both are mobile + can be multi-modal
- Spatial resolution → EEG 5-9cm, fNIRS 1-3cm
- Temporal resolution → EEG >1000Hz, fNIRS 10Hz
- EEG can be time locked to each step
- fNIRS can be time locked to blocks of steps
- EEG is sensitive to electrical noise in environment
- Both modalities susceptible to movement artifacts → fNIRS may have better resistance to movement artifact
- ***Which modality to choose? Depends on your research question***

Example Experiment

- Aim → use fNIRS to delineate brain activation differences between 'Active' and 'Passive' overground gait in a robotic exoskeleton
- Multimodal → fNIRS + EMG + Xsens + robotic exoskeleton (Ekso)




Sharable link to article: <https://rdcu.be/esed3>

(Peters et al. *J Neuroeng & Rehab* 2020)


Example Experiment

- Montage → coverage of brain regions that relate to your research question. In Ekso study, wanted broad coverage of frontal and parietal regions bilaterally



Shareable link to article: <https://rdcu.be/esed3>
(Peters et al. *J Neuroeng & Rehab* 2020)






sue.peters@uwo.ca

Example Experiment

- Montage → include short channels



Review (Chen et al. Front Neurosci 2020)

Sharable link to article: <https://rdcu.be/esed3>

(Peters et al. *J Neuroeng & Rehab* 2020)



Example Experiment

- Multimodal synchronization with fNIRS
 - EMG, IMU (Xsens)
 - Triggers during data collection across modalities are required → LSL or TTL
- Could adapt this multimodal synchronization
 - We will try fNIRS with WINGS today
 - Synchronize fNIRS with ECG, EMG, GSR, etc

Sharable link to article: <https://rdcu.be/esed3>

(Peters et al. *J Neuroeng & Rehab* 2020)



Example Experiment

- Set up tips (i.e., didn't make it into the publication....)
 - fNIRS cap placed AFTER the EMG and IMUs place and measurements completed for fitting the robot
 - Used different gromet tops with low tension ('zero'-low tension springs) on the forehead
 - fNIRS unit worn by trainee triggering the robot
 - Check +++ to ensure the cap/wires secure
- Motion artifact: Tried Wavelet and Spline and Wavelet+Spline....chose Wavelet
- Consecutive gaussian function for HRF estimation

Sharable link to article: <https://rdcu.be/esed3>

(Peters et al. *J Neuroeng & Rehab* 2020)

Methods for patient populations


Post-stroke gait

- Reduced velocity
- Decreased step length, coordination
- Is fatiguing → impacts to study design
- Channels grouped into brain regions
 - Option to use 3D digitization
- Need all lesions on the same side for statistics
 - Some brains 'flipped'
 - Results: ipsilesional & contralesional activation (instead of R and L)



Study design – depends on research question

- Design
 - Block (repeating the same task) vs. Event (task done once)



ST = single task, DT = dual task

Study design – depends on research question

- Analysis methods

- GLM

- Makes assumptions about the shape of the HRF
 - Common to assume the shape is fixed
 - Canonical, boxcar, gaussian (consecutive sequence)


- Correlations

- Between brain regions
 - Common in resting state

- Deconvolution

- Allows the shape of the HRF to vary
 - Pros: Allows participant/brain region variability
 - Cons: need lots of data, freedom in HRF shape makes interpretation more difficult (may be okay if underlying neural and hemodynamics of population well understood)


- Other?



Study design - considerations for balance and gait research questions


Location of leg region in M1

- Motor homunculus
 - Leg, foot, etc located in the midline/longitudinal fissure
- Gait is a whole-body task
 - postural control, arm swing
 - +/- cognitive aspects of gait
- Consider a control condition in standing
 - Examples
 - Quiet standing
 - Standing with a cognitive task
 - Standing with button press



Tips & Tricks for using fNIRS during Gait

- Use timing and live observation of signal to determine when to start gait trial
- Gait
 - Long hallway (measure the distance), Treadmill, Outdoors?
- Quiet environment
 - (unless studying auditory or cognitive effects of distraction?)
- Cannot clean data that is very noisy → physical set up is key!!
 - Observe heart rate
 - Wires that 'bounce' during gait generate noise
 - Use multiple holders + Velcro + tighten
- Short channels
 - Use multiple holders to secure so these wires don't overlap with holes for optodes



Hands on


Wrap up: fNIRS tips for setting up multimodal research projects

- Study design

- Consider your participants → what can they tolerate? How much data can you reasonably collect in a session?
 - Number of blocks, signal to noise, etc
- Patients → what other measures do you want to consider (kinematics? physiology including ECG, respiration rate etc, EMG?)

- Analysis

- Measure a time window around the HRF peak? Or the whole block?
- GLM vs correlation vs other?



(Lim et al. Front Neurol 2022) 20

Wrap up: fNIRS tips for setting up multimodal research projects

- Collaborators can help – consider researchers and industry
- Lots of piloting of data collection and analysis
- Start with replication - gait vs rest, finger tapping vs rest
- Read recommendation publications → these are not specific to multimodal data collection but a good place to start
 - Menant et al. 2020: A consensus guide to using functional near-infrared spectroscopy in posture and gait research
 - Yücel et al. 2021: Best practices for fNIRS publications
- Field is evolving very quickly!
- Expect more recommendations/guideline papers and updates to current ones

Wrap up: EEG tips for setting up multimodal research projects

- Many similarities with fNIRS tips
- Consider ways to minimize noise, use noise reference channels
- Consider what system to use
 - If your own system isn't possible – online data sets
 - <https://openneuro.org>
 - <https://nemar.org>
- Include MOBI in your research!

Thank you

Questions?

Website with all slides +
additional materials

